Synthes ProPlan CMF

Deutsch

English

Español

Français

Italiano

Português

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Acrylic ProPlan CMF Orthognathic Splints

This document contains general instructions for use for ProPlan CMF Orthognathic Splints. For case-specific instructions, refer to the case report.

DESCRIPTION

ProPlan CMF Orthognathic Splints are custom-made devices designed to fit, or represent, the patient's anatomy. They are intended for improving and simplifying the performance of surgical interventions.

INDICATIONS

ProPlan CMF Orthognathic Splints are intended to be used as surgical tools to transfer a pre-operative plan to surgery. ProPlan CMF Orthognathic Splints are intended to guide the alignment of the maxilla and / or mandible in orthognathic surgical procedures. ProPlan CMF Orthognathic Splints are intended for single use only.

MATERIAL

Hybrid epoxy-acrylate polymer

CONTRAINDICATIONS

Do not use in the case of active infection of the surgical area.

STORAGE

It is advised to store the splints at room temperature in a properly cleaned, dry and dark place. Avoid exposing the splints to UV-light. Only open the package right before preparing the splints for surgery (i.e. before cleaning and sterilization).

WARNING

- The user should be aware of possible allergic reactions to materials used in the splints. The patient should be informed on this matter by the user.
- These are patient specific, single use, disposable splints.
- Do not attempt to reuse or recondition the splints.
- Do not alter the splints in any way.
- The splints are to be used by a trained physician in the performance of surgery.
- Be aware that these patient specific splints have been manufactured based on CT/MRI scans of the patient and plaster models of the patient's teeth. If the patient's anatomy has changed significantly since the CT/MRI scan or the moment patient's tooth impressions were taken, the splints should not be used. This includes teeth movement due to ongoing orthodontic treatment.
- The splints should be properly cleaned before sterilization. Do not use if they are broken, cracked, or are visibly contaminated.
- The splints in this package are provided non-sterile.

Distributed by: Synthes GmbH
Eimattstrasse 3
4436 Oberdorf
Switzerland

Manufactured by: Materialise NV
Technologielaan 15
3001 Leuven
Belgium

PRECAUTIONS

- The splints are intended to be used during the surgery only. If, for any reason, it is decided to keep a splint in the patient's mouth for a longer period, this should never be longer than 24 hours.
- It is advised to use the splints within 6 months of performing the CT/MRI scans on which they are based. If the patient's anatomy has changed significantly since the time of the CT/MRI-scan, the splints should not be used, even if the time period of 6 months has not expired.
- Do not apply excessive force on the splints, or place heavy objects on top.
- Markings on the splints used to indicate anatomical references and case information must be legible. These include lines indicating anatomical directions, identifiers with case information such as implant size, and the unique case identifier (see below). Notify your representative if the markings are not legible or if the identifiers do not correspond to the intended patient or surgeon.

PATIENT SPECIFIC SPLINT IDENTIFIERS

A unique identifier is indicated on each splint. This alphanumeric code corresponds to the last 6 digits of the case identifier. The case report shipped with each patient case details the link between the case identifier and the patient's identity and the shipped instruments.

Before using the splint, check the unique identifier for readability and confirm that it corresponds with the case identifier.

POSSIBLE ADVERSE EFFECTS

Infection following the surgical procedure. Introduction of foreign materials can result in an inflammatory response or allergic reaction.

INSTRUCTIONS FOR USE

• Fitting of the splint

- The splint is designed to fit the patients' anatomy. The supporting surface (teeth, mucosa) should be completely freed to assure proper fit of the splint.
- Limit the forces on the splint to 20N (i.e. 4.5 lb(f)).
- Take enough time to fit the splint on the patient. Try different positions and check whether or not the splint stays in place. Choose the most stable position, i.e. the position in which the least movement of the splint is possible. Don't push the splint down too hard to avoid breaking the splint.

If it is not possible to place the splint on the patient in a unique and stable position, the splint does not guarantee an accurate transfer of the pre-operative planning.

- Do not alter the splint before use. Small particles might come off, which could contaminate the operating region. In addition, altering the size of the splint may result in an inadequate fit to the patient's anatomy or reduced mechanical strength of the splint. Therefore, it is the responsibility of the user if the splint is altered in any way prior to, or during, surgery.

• During splint use

- Make sure the splint maintains its position between maxilla and mandible during surgery.
- In case a sequence of multiple splints is provided, start with the splint whose label ends in the number '1'. Work your way up from there, using splints with increasing label numbers until you've used all the splints.

CLEANING AND STERILIZATION INSTRUCTIONS

ProPlan CMF Splints are NOT STERILE and must be thoroughly cleaned and sterilized prior to use in surgery

CLEANING

Whenever possible, a washer/disinfector (according to ISO 15883) and ultrasonic cleaning equipment should be used to clean the splints. The detergents and/or enzymatic cleaner should have a pH value in the range pH 7-11. The splints can be cleaned using manual cleaning and/or automated cleaning in a washer/disinfector with manual pre-cleaning and ultrasonic cleaning.

Manual cleaning:

Step	Cleaning instructions
1	Prepare a fresh, newly-made solution using warm de-ionized (DI) or purified water (PURW) and enzymatic cleaner or detergent.
2	Carefully wash the splint manually.
3	Rinse the splint thoroughly with DI or PURW.
4	Dry the splint using a clean, soft, lint-free cloth or clean compressed air.

Manual pre-cleaning:

Step	Minimum Duration	Cleaning instructions
1	1 minute	Rinse the splint under running cold tap water.
2	2 minutes	Manually clean the splint in a newly-made enzymatic cleaner or detergent solution.
3	1 minute	Rinse the splint using cool to lukewarm running tap water. Use a syringe, pipette or water pistol to flush cylinders, slots, and other hard-to-reach areas.
4	15 minutes	Clean the splint ultrasonically per manufacturer's recommended temperature (usually 32°-60°C or 90°-140°F) and specially formulated detergents. Follow manufacturer's recommendations for proper cleaning solution formulated specifically for ultrasonic cleaners and medical equipment.
5	2 minutes	Rinse the splint using DI or PURW. Use a syringe, pipette, or water pistol to flush cylinders, slots, and other hard-to-reach areas.

Automated cleaning in a washer/disinfector:

Step	Minimum Duration	Cleaning instructions
Pre-wash	2 minutes	Cold tap water
Wash	10 minutes	Warm tap water (>40°C); use detergent
Neutralize	2 minutes	Warm tap water with neutralizer, if necessary
Rinse	2 minutes	Rinse with warm DI or PURW (>40°C)
Thermal disinfection	7 minutes	At minimum 94°C
Dry	40 minutes	At minimum 90°C

Before the cleaned products are packaged and sterilized, carefully examine them to see if they are clean and undamaged.

STERILIZATION

Recommended sterilization specifications.

The splints can be steam sterilized up to two (2) times prior to use. Users should conduct testing in the health care facility to ensure that conditions essential to sterilization can be achieved.

Sterilize the splints using pre-vacuum steam sterilization before use.

During sterilization of single models pouches may be used. Only standard medical grade steam sterilization polyethylene or Tyvek pouches should be used. Ensure that the pouch is large enough to contain the splint without stressing the seals or tearing the pouch.

Don't exert force on the splints right after autoclaving, as the splints might deform during the cool-down period. After sterilization, the color of the splints can change (yellow/brown), but this does not affect their properties.

Use one of the following standard steam sterilization settings:

- Pre-vacuum Cycle Worldwide^{2,3}:
 - Minimum temperature: 132°C
 - Maximum temperature: 135°C
 - Minimum exposure time: 4 minutes
 - Maximum exposure time: 5 minutes
 - Minimum vacuum drying time: 30 minutes

- Pre-vacuum Cycle UK, NL^{1,3}:
 - Minimum temperature: 134°C
 - Maximum temperature: 135°C
 - Minimum exposure time: 3 minutes
 - Maximum exposure time: 5 minutes
 - Minimum vacuum drying time: 30 minutes

1 Minimum validated steam sterilization time required to achieve a 10⁻⁶ sterility assurance level (SAL).

2 Minimum validated steam sterilization temperature required to achieve a 10⁻⁶ sterility assurance level (SAL).

3 In the case local or national specifications for steam sterilization requirements are stricter or more conservative than those listed in this table, please contact Materialise before sterilizing and using the splints.

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